

Product Datasheet

PRUNE2 antibody (orb19884)

Description	Goat polyclonal antibody to PRUNE2
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PRUNE2 / BMCC1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 30kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 29.5kDa according to NP_620173.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	MESEKISEKQEE
Clonality	Polyclonal
MW	29.5
Entrez	158471
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19884
(1µg/ml)
staining of
Human Brai...