

## Product Datasheet

# Anti-Histone H3 (tri methyl K27) Rabbit Monoclonal Antibody (orb1993265)

<b>Description</b>	Anti-Histone H3 (tri methyl K27) Rabbit Monoclonal Antibody catalog # orb1993265. Tested in WB, IHC,...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC, WB
<b>Immunogen</b>	A synthesized peptide derived from human Histone H3 (tri methyl K27)
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	WB 1:500-1:2000IHC 1:50- 1:200ICC/IF 1:50-1:200
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P68431</b>
<b>Expiration Date</b>	12 months from date of receipt.