

Product Datasheet

Anti-Histone H3 (di methyl K79) Rabbit Monoclonal Antibody (orb1993283)

Description	Anti-Histone H3 (di methyl K79) Rabbit Monoclonal Antibody catalog # orb1993283. Tested in WB, IHC, ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	A synthesized peptide derived from human Histone H3 (di methyl K79)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IHC 1:50- 1:200ICC/IF 1:50-1:200
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.