

Product Datasheet

Anti-HRPT2 / CDC73 Rabbit Monoclonal Antibody (orb1993313)

Description	Anti-HRPT2 / CDC73 Rabbit Monoclonal Antibody catalog # orb1993313. Tested in WB applications. This ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human HRPT2 / CDC73
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q6P1J9
Expiration Date	12 months from date of receipt.