

Product Datasheet

CD123 Antibody (PerCP/Cy5.5) [Out of stock] (orb1993392)

Description	CD123 is the alpha chain of interleukin 3 receptor (IL-3R alpha). This subunit heterodimerizes with ...
Reactivity	Human
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	IL3 receptor alpha chain expressed on the surface of transiently transfected COS cells
Target	CD123
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P26951
Entrez	3563
Expiration Date	12 months from date of receipt.