

Product Datasheet

CD16 Antibody (PerCP/Cy5.5) (orb1993405)

Description	CD16 (FcγRIII) is a 50-65 kDa glycoprotein serving as a low affinity IgG receptor. Human Fcγma...
Reactivity	Human, Primate
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	Human neutrophils
Target	CD16
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.