

Product Datasheet

CD34 Antibody (PerCP/Cy5.5) (orb1993419)

Description	CD34 is a highly glycosylated monomeric 111-115 kDa surface protein, which is present on many stem c...
Reactivity	Human, Primate
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	Human endothelial vesicles
Target	CD34
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P28906
Entrez	947
Expiration Date	12 months from date of receipt.