

## Product Datasheet

### CD4 Antibody (PerCP/Cy5.5) (orb1993427)

<b>Description</b>	CD4 (T4) is a single chain transmembrane glycoprotein and belongs to immunoglobulin supergene family...
<b>Reactivity</b>	Human
<b>Conjugation</b>	PerCP/Cy5.5
<b>Tested Applications</b>	FC
<b>Immunogen</b>	2 N-terminal domains of human CD4 fused to human IgG1 Fc
<b>Target</b>	CD4
<b>Preservatives</b>	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
<b>Storage</b>	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
<b>Note</b>	For research use only
<b>Application notes</b>	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P01730</a>
<b>Entrez</b>	<a href="#">920</a>
<b>Expiration Date</b>	12 months from date of receipt.