

Product Datasheet

CD25 Antibody (PerCP/Cy5.5) (orb1993438)

Description	CD25 (IL2Ralpha, Tac) is a ligand-binding alpha subunit of interleukin 2 receptor (IL2R). Together w...
Reactivity	Human
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	PHA-activated peripheral blood leucocytes
Target	CD25
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P01589
Entrez	3559
Expiration Date	12 months from date of receipt.