

Product Datasheet

CD152 Antibody (PE/Cy7) (orb1993503)

Description	CD152 / CTLA-4 is a homodimeric transmembrane protein similar to CD28 and binding the same ligands, ...
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	Human CD152-IgG heavy chain fusion protein
Target	CD152
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Uniprot ID	P16410
Entrez	1493
Expiration Date	12 months from date of receipt.