

Product Datasheet

EGFR Antibody (PE/Cy7) (orb1993524)

Description	The oncoprotein EGFR (epidermal growth factor receptor), also known as HER1 / ErbB1, is a member of ...
Reactivity	Equine, Human
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	Human epidermoid carcinoma line A431
Target	EGFR
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Uniprot ID	P00533
Entrez	1956
Expiration Date	12 months from date of receipt.