

Product Datasheet

Ig Kappa Light Chain Antibody (PE/Cy7) (orb1993527)

Description	Immunoglobulin classes share the same basic four polypeptide chain structure of two heavy chains (fi...
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	Human IgG-kappa myeloma protein
Target	Ig Kappa Light Chain
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests. Extracellular and intracellular staining.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.