

Product Datasheet

CD28 Antibody (PE/Cy7) (orb1993548)

Description	CD28 is the critical T cell costimulatory receptor which provides to the cell the important second a...
Reactivity	Human, Primate
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	DC28.1.3.3 murine T cell hybridoma transfected with human CD28 cDNA
Target	CD28
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	P10747
Entrez	940
Expiration Date	12 months from date of receipt.