

Product Datasheet

CD86 Antibody (PE/Cy7) (orb1993551)

Description	CD80 (B7-1) and CD86 (B7-2) are ligands of T cell critical costimulatory molecule CD28 and of an inh...
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	B-lymphoblastoid cell line ARH 77
Target	CD86
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P42081
Entrez	942
Expiration Date	12 months from date of receipt.