

Product Datasheet

CD273 Antibody (PE/Cy7) (orb1993579)

Description	CD273 / PD-L2 (programmed death ligand-1), also known as B7-DC, is a member of the B7 family of regu...
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	human CD273 cDNA
Target	CD273
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Uniprot ID	Q9BQ51
Entrez	80380
Expiration Date	12 months from date of receipt.