

Product Datasheet

CD279 Antibody (PE/Cy7) (orb1993581)

Description	CD279 / PD-1 (programmed cell death 1), a transmembrane protein of CD28/CTLA-4 family. It is express...
Reactivity	Human, Primate
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	recombinant human CD279
Target	CD279
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	Q15116
Entrez	5133
Expiration Date	12 months from date of receipt.