

Product Datasheet

CD49d Antibody (APC/Cy7) (orb1993630)

Description	CD49d / integrin alpha 4, unlike other alpha integrins, neither contains an I-domain, nor undergoes ...
Reactivity	Bovine, Canine, Equine, Feline, Human, Primate, Sheep
Conjugation	APC/Cy7
Tested Applications	FC
Target	CD49d
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	P13612
Entrez	3676
Expiration Date	12 months from date of receipt.