

Product Datasheet

HLA-DR Antibody (Pacific Orange) [Out of stock] (orb1993635)

Description	HLA-DR, a member of MHC class II glycoproteins, that bind intracellularly processed peptides and pre...
Reactivity	Canine, Human, Primate
Conjugation	Pacific Orange
Tested Applications	FC
Immunogen	Human B lymphocytes
Target	HLA-DR
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.