



Product Datasheet

CD34 Antibody (PerCP) (orb1993652)

Description	CD34 is a highly glycosylated monomeric 111-115 kDa surface protein, which is present on many stem c...
Reactivity	Human
Conjugation	PerCP
Tested Applications	FC
Immunogen	Permanent human cell line derived from peripheral leucocytes of a patient suffering from chronic myeloid leukaemia.
Target	CD34
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P28906
Entrez	947
Expiration Date	12 months from date of receipt.



Separation of human CD34 positive stem c...



Flow cytometry surface staining pattern ...