

Product Datasheet

CD235a Antibody (Pacific Blue) (orb1993660)

Description	CD235a (glycophorin A, GPA) is a transmembrane sialoglycoprotein expressed on erythrocytes and their...
Reactivity	Human, Rat
Conjugation	Pacific Blue
Tested Applications	FC
Immunogen	Membrane preparation from splenic hairy cell leukemia
Target	CD235a
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P02724
Entrez	2993
Expiration Date	12 months from date of receipt.