

Product Datasheet

CD300e Antibody (Pacific Blue) (orb1993662)

Description	CD300e / IREM-2 (immune receptor expressed by myeloid cells 2), also known as CLM2 or LMIR6, is a mo...
Reactivity	Human
Conjugation	Pacific Blue
Tested Applications	FC
Immunogen	CD300e-HA-transfected cells
Target	CD300e
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	Q496F6
Entrez	342510
Expiration Date	12 months from date of receipt.