

Product Datasheet

CD6 Antibody (Pacific Blue) (orb1993724)

Description	CD6, also known as T12, is a member of the scavenger receptor superfamily found on T and B cell subs...
Reactivity	Human
Conjugation	Pacific Blue
Tested Applications	FC
Immunogen	Human CD6 antigen purified by immunoaffinity chromatography from HBP-ALL cells followed by preparative SDS-PAGE of non-boiled non-reduced sample (excised piece of gel corresponding to the 100 kDa zone).
Target	CD6
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P30203
Entrez	923
Expiration Date	12 months from date of receipt.