



Product Datasheet

CD97 Antibody (PE) (orb1993766)

Description	CD97 is a G-protein-coupled seven-span transmembrane adhesive receptor that is constitutively expres...
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	PHA-activated peripheral blood cells
Target	CD97
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P48960
Entrez	976
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-fille...



Flow cytometry surface staining pattern ...