

Product Datasheet

CD1a Antibody (PE) (orb1993768)

Description	CD1a, together with CD1b and c, belongs to group 1 of CD1 glycoproteins. These proteins serve as ant...
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	Human thymocytes
Target	CD1a
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Uniprot ID	P06126
Entrez	909
Expiration Date	12 months from date of receipt.