

Product Datasheet

CD1a Antibody (APC) (orb1993807)

| | |
|----------------------------|---|
| Description | CD1a, together with CD1b and c, belongs to group 1 of CD1 glycoproteins. These proteins serve as ant... |
| Reactivity | Human |
| Conjugation | APC |
| Tested Applications | FC |
| Immunogen | Human thymocytes |
| Target | CD1a |
| Preservatives | Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide |
| Storage | Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze. |
| Note | For research use only |
| Application notes | Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. |
| Isotype | Mouse IgG2b kappa |
| Clonality | Monoclonal |
| Uniprot ID | P06126 |
| Entrez | 909 |
| Expiration Date | 12 months from date of receipt. |