

## Product Datasheet

# Rabbit anti-Goat IgG (H&L) - Affinity Pure, HRP Conjugate, Cross-Adsorbed (orb1993898)

<b>Description</b>	Rabbit anti-Goat IgG (H&L) - Affinity Pure, HRP Conjugate, Cross-Adsorbed
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Goat
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.