

Product Datasheet

Donkey anti-Goat IgG (H&L) - Affinity Pure, HRP Conjugate, Cross-Adsorbed (orb1993915)

| | |
|----------------------------|---|
| Description | Donkey anti-Goat IgG (H&L) - Affinity Pure, HRP Conjugate, Cross-Adsorbed |
| Species/Host | Donkey |
| Reactivity | Goat |
| Conjugation | HRP |
| Tested Applications | ELISA, IHC, WB |
| Note | For research use only |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |