

Product Datasheet

Goat anti-Rabbit IgG (H&L) - Affinity Pure, FITC Conjugate, Cross-Adsorbed (orb1993966)

| | |
|----------------------------|--|
| Description | Goat anti-Rabbit IgG (H&L) - Affinity Pure, FITC Conjugate, Cross-Adsorbed |
| Species/Host | Goat |
| Reactivity | Rabbit |
| Conjugation | FITC |
| Tested Applications | ELISA, FC, IF, IM |
| Note | For research use only |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |