

## Product Datasheet

# violetFluor 500 Mouse IgG2a Isotype Control (C1.18.4) Antibody (orb1994120)

<b>Description</b>	violetFluor 500 Mouse IgG2a Isotype Control (C1.18.4) Antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	violetFluor 500
<b>Tested Applications</b>	FC
<b>Storage</b>	2-8°C protected from light
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG2a
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.