

## Product Datasheet

### PE Anti-Human CD200 (OX-104) (orb1994391)

<b>Description</b>	PE Anti-Human CD200 (OX-104)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	PE
<b>Tested Applications</b>	FC
<b>Storage</b>	2-8°C protected from light
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG1, kappa
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P41217</a>
<b>Entrez</b>	<a href="#">4345</a>
<b>Expiration Date</b>	12 months from date of receipt.