


Product Datasheet

IRAK3 antibody (orb19959)

Description	Goat polyclonal antibody to IRAK3
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IRAK3 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx. 65kDa band observed in Rat Liver lysates (calculated MW of 67.8kDa according to NP_082955.2). Recommended concentration: 0.05-0.2 µg/ml.
Protein Sequence	GHSYGSKPMEKR
Clonality	Polyclonal
MW	67.8
Dilution Range	ELISA: 1:64000, WB: 0.05-0.2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Liver lysat...



Western blot analysis of Mouse Bone Marr...