



## Product Datasheet

### UXT antibody (orb19973)

<b>Description</b>	Goat polyclonal to UXT
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Target</b>	UXT
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody is expected to recognize both reported isoforms (NP_004173.1; NP_705582.1).ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 20kDa band observed in Human Kidney and Skeletal Muscle lysates (calculated MW of 19.6kDa according to NP_705582.1). Recommended concentration: 0.01-0.03 µg/ml.
<b>Protein Sequence</b>	DHRDKVYEQLAKY
<b>Clonality</b>	Polyclonal
<b>MW</b>	18.2; 19.6
<b>Entrez</b>	<b>8409</b>
<b>Dilution Range</b>	ELISA: 1:64000, WB: 0.01-0.03 µg/ml, IHC-P: 5-10µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Human Kidney ly...



Immunohistochemical staining of Human Th...