

## Product Datasheet

### Slc10a2\* antibody (orb19974)

<b>Description</b>	Goat polyclonal antibody to Slc10a2*
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	Slc10a2 (mouse)
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 35kDa band observed in Mouse Small Intestine lysates (calculated MW of 38.1kDa according to NP_035518.1). Recommended concentration: 0.5-1.5 µg/ml.
<b>Protein Sequence</b>	DETNKGFPDEK
<b>Clonality</b>	Polyclonal
<b>MW</b>	38.1
<b>Dilution Range</b>	ELISA: 1:32000, WB: 0.5-1.5 µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



orb19974  
(0.5ug/ml)  
staining of  
Mouse Sm...