

Product Datasheet

UBE2A Peptide - N-terminal region (orb1997751)

Description	UBE2A Peptide - N-terminal region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DFKRLQEDPPAGVSGAPSENNIMVWNAVIFGPEGTPFEDGTFKLTIEFTE
MW	17 kDa
Uniprot ID	Q9Z255
NCBI	NP_001269090.1
Expiration Date	6 months from date of receipt.