

Product Datasheet

USP25 Peptide - N-terminal region (orb1998321)

Description	USP25 Peptide - N-terminal region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: AISRVLEASIAENKACLKRTPEVWRDSRNPYDRKRQDKAPVGLKNVGNT
MW	123 kDa
Uniprot ID	Q9UHP3
NCBI	NP_001269970.1
Expiration Date	6 months from date of receipt.