

Product Datasheet

UBA7 Peptide - middle region (orb1998852)

Description	UBA7 Peptide - middle region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: QEVHHAHCLHQAFKALHKFQHLHGRPPQPWDPVDAETVVGLARDLEPLKR
MW	111 kDa
Uniprot ID	P41226
NCBI	NP_003326.2
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)