

Product Datasheet

PRDM16 Peptide - middle region (orb1999239)

Description	PRDM16 Peptide - middle region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SAPASGEEQPLDLSIGSRARASQNGGGREPRKNHVYGERKLGAGEGLPQV
MW	120 kDa
Uniprot ID	Q9HAZ2
NCBI	NP_071397.3
Expiration Date	6 months from date of receipt.