

Product Datasheet

UBE2E3 Peptide - N-terminal region (orb1999424)

| | |
|-------------------------|---|
| Description | UBE2E3 Peptide - N-terminal region |
| Preservatives | Lyophilized powder |
| Form/Appearance | Lyophilized powder |
| Storage | Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles. |
| Note | For research use only |
| Protein Sequence | Synthetic peptide located within the following region: MSSDRQRSDDESPSTSSGSSDADQRDPAAPEPEEQEERKPSATQQKKNTK |
| MW | 49 kDa |
| Uniprot ID | Q8WV16 |
| NCBI | NP_001265483.1 |
| Expiration Date | 6 months from date of receipt. |