

Product Datasheet

MCCC2 Peptide - C-terminal region (orb1999792)

Description	MCCC2 Peptide - C-terminal region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: QFSSADEAALKEPIIKKFE EEGNPYYSSARVWDDGIIDPADTRLVLGLSF
MW	57 kDa
Uniprot ID	Q9HCC0
NCBI	NP_071415.1
Expiration Date	6 months from date of receipt.