

Product Datasheet

ITIH2 Peptide - C-terminal region (orb1999835)

Description	ITIH2 Peptide - C-terminal region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: NVDFLGIYIPPTNKFSPKAHGLIGQFMQEPKIHIFNERPGKDPEKPEASM
MW	106 kDa
Uniprot ID	P19823
NCBI	NP_002207.2
Expiration Date	6 months from date of receipt.