

Product Datasheet

NVL Peptide - N-terminal region (orb1999889)

Description	NVL Peptide - N-terminal region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	84 kDa
Uniprot ID	O15381
NCBI	NP_001230075.1
Expiration Date	6 months from date of receipt.