

Product Datasheet

MMACHC Peptide - C-terminal region (orb2000080)

Description	MMACHC Peptide - C-terminal region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KAYFSTPPAQLALLGLAQPSEKSSPDLPFTTPAPKKPGNPSRARSW
MW	31 kDa
Uniprot ID	Q9Y4U1
NCBI	NP_056321.2
Expiration Date	6 months from date of receipt.