

## Product Datasheet

### USP19 Peptide - middle region (orb2000463)

<b>Description</b>	USP19 Peptide - middle region
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	This is a synthetic peptide designed for use in combination with USP19 Rabbit Polyclonal Antibody (orb589350). It may block above mentioned antibody from binding to its target protein in western blot and/or immunohistochemistry under proper experimental settings.
<b>Protein Sequence</b>	Synthetic peptide located within the following region: SQFTGYAQHDAQEFMAFLLDGLHEDLNRIQNKPYTETVDS DGRPDEVVAE
<b>MW</b>	156 kDa
<b>Uniprot ID</b>	<b>O94966</b>
<b>NCBI</b>	<b>NP_001186089.1</b>
<b>Expiration Date</b>	6 months from date of receipt.

---

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)