

Product Datasheet

PXMP2 Peptide - middle region (orb2000623)

Description	PXMP2 Peptide - middle region
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	21 kDa
Uniprot ID	Q9NR77
NCBI	NP_061133.1
Expiration Date	6 months from date of receipt.