

## Product Datasheet

# APOBEC3C Peptide - C-terminal region (orb2000646)

<b>Description</b>	APOBEC3C Peptide - C-terminal region
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Synthetic peptide located within the following region: EGLRSLSEQGVAVEIMDYEDFKYCWENFVYNDNEPFKPKWGLKTNFRLK
<b>MW</b>	23 kDa
<b>Uniprot ID</b>	<a href="#">Q9NRW3</a>
<b>NCBI</b>	<a href="#">NP_055323.2</a>
<b>Expiration Date</b>	6 months from date of receipt.