

Product Datasheet

CNTN2 Peptide - middle region (orb2000895)

Description	CNTN2 Peptide - middle region
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LGDNLTLQCHASHDPTMDLTFTWTLDDFPIDFDKPGGHYRRTNVKETIGD
MW	114 kDa
Uniprot ID	Q02246
NCBI	NP_005067.1
Expiration Date	6 months from date of receipt.