

Product Datasheet

APOBEC3C Peptide - C-terminal region (orb2001233)

Description	APOBEC3C Peptide - C-terminal region
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EGLRSLSQEGVAVEIMDYEDFKYCWENFVYNDNEPFKPKWGLKTNFRLK
MW	23 kDa
Uniprot ID	Q9NRW3
NCBI	NP_055323.2
Expiration Date	6 months from date of receipt.