

## Product Datasheet

### TNIK Peptide - middle region (orb2001431)

<b>Description</b>	TNIK Peptide - middle region
<b>Tested Applications</b>	WB
<b>Preservatives</b>	Lyophilized powder
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Synthetic peptide located within the following region: ETSHADSFSGSISREGTLMIRETSGEKKRSGHSDSNGFAGHINLPDLVQQ
<b>MW</b>	155 kDa
<b>Uniprot ID</b>	<a href="#">Q9UKE5</a>
<b>NCBI</b>	<a href="#">NP_001155032.1</a>
<b>Expiration Date</b>	6 months from date of receipt.