

Product Datasheet

LGI3 Peptide (orb2001688)

Description	This is a synthetic peptide designed for use in combination with anti-LGI3 antibody (Catalog #: ARP5...
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	60
Uniprot ID	Q8N145
NCBI	NP_644807.1
Expiration Date	6 months from date of receipt.