

Product Datasheet

EHHADH Peptide - middle region (orb2001747)

Description	This is a synthetic peptide designed for use in combination with anti-ECHP Antibody (ARP76113_P050),...
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VIAVDSDKNQLATANKMITSVLEKEASKMQSGHPWSGPKPRLTSSVKEL
MW	79kDa
Uniprot ID	Q08426
NCBI	NP_001957
Expiration Date	6 months from date of receipt.